

OPHTHALMIC SURGICAL INSTRUMENTS



Version 250108

Contents

Premium Instruments.....	3
Forceps.....	3
Tying Forceps.....	11
Toothed Forceps.....	12
Chalazion Forceps.....	17
Meibomian Forceps.....	18
Lens Forceps.....	19
Conjunctiva Forceps.....	20
Prechopper.....	20
Capsulorhexis Forceps.....	22
Miscellaneous Forceps.....	26
Scissors.....	19
Lens Cutter.....	21
Iris Scissors.....	18
Corneal Scissors.....	29
Tenotomy & Utility Scissors.....	32
Vannas Scissors.....	33
Suture-Stitch Scissors.....	36
Conjunctiva Scissors.....	36
Needle Holders.....	37
Hooks Spatulas & Manipulators.....	40
Strabismus Hooks & Retractors.....	45
Phaco Choppers.....	47
Speculum.....	52
ICL ® Instruments.....	56
LASIK LASEK & ReLEx SMILE.....	57
Markers.....	61
Irrigation & Aspiration.....	64
Calipers.....	67
Punch.....	68
Lacrimal Instruments & Trabeculotome Set.....	68
Cannulas.....	69
Corneal Ring Set.....	72
Mermoud Set.....	73
MISCELLANEOUS.....	74
Injectors & Corneal Trephine.....	78
Sapphire Knives.....	79
Sterilizing Trays.....	81
Vitreo-Retinal 20G/23G Instruments.....	82
Vitreo-Retinal 25G Instruments.....	85
Micro Capsulorhexis Forceps.....	86
Micro IOL Removal Forceps & Scissors.....	87




Premium Forceps

McPherson Tying Forceps

Premium design, 5mm tying platforms
82mm total length


PM1181 Straight Shaft
PM1182 Curved Shaft
PM1183 Angled Shaft



McPherson Tying Forceps

Premium design, 5mm tying platforms
105mm total length

PM1184 Straight Shaft
PM1185 Curved Shaft
PM1186 Angled Shaft



Bonn Toothed Forceps

Premium design, 5mm tying platforms
0.12mm 1 × 2 teeth, 82mm total length

PM1442 Straight Shaft
PM1443 Curved Shaft
PM1444 Angled Shaft



Premium Forceps

Bonn Toothed Forceps

Premium design, 5mm tying platforms
0.12mm 1×2 teeth, 105mm total length

- PM1445 Straight Shaft
- PM1446 Curved Shaft
- PM1447 Angled Shaft



Bonn Toothed Forceps

Premium design, 5mm tying platforms
Straight Shaft, 105mm total length

- PM1448 0.3mm 1×2 teeth
- PM1449 0.4mm 1×2 teeth
- PM1450 0.5mm 1×2 teeth



Tying Forceps

Premium design, 5mm tying platforms, 110mm total length
Made of Ultra high hardness Stainless Steel

- PM1452 Straight Shaft



Toothed Forceps

Premium design, 5mm tying platforms, 0.12mm 1×2 teeth, 110mm total length
Made of Ultra high hardness Stainless Steel

- PM1453 Straight Shaft



Premium Scissors

Westcott Curved Scissors

Premium design, Made of Ultra high hardness Stainless Steel

TN3332 Sharp tips, 18mm from pivot to tip, 130mm length

TN3333 Blunt tips, 18mm from pivot to tip, 130mm length



Premium Scissors

Vannas Curved Capsulotomy Scissors

Premium design, Made of Ultra high hardness Stainless Steel

VN3444 Sharp pointed tips 16mm from pivot to tip, 110mm length



Vannas Straight Scissors

Made of Ultra high hardness Stainless Steel

PM3423 Sharp pointed tips, 8.5mm blades,
85mm total length



Vannas Curved Scissors

Made of Ultra high hardness Stainless Steel

PM3424 Sharp pointed tips, 8.5mm blades,
85mm total length



Vannas Angled Scissors

Made of Ultra high hardness Stainless Steel

PM3425 Sharp tips, 6.5mm blades,
85mm total length



Premium Phaco Tip

Phaco Tip & Wrench

PM2331 Wrench for Phaco Tip



Premium Needle Holders

Castroviejo Straight Needle Holder

Premium design, Cross action shaft

NH4301 Cross action 10mm jaws, blunt, without lock, 120mm length



Premium Phaco Choppers

Nagahara Phaco Chopper

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM6113 Angled Curved shaft, 1.1mm tip

PM6114 Angled Curved shaft, 1.6mm tip



Seibel Nucleus Chopper

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM6192 Chopper with ball-shaped tip for efficient capsule protection



Premium Phaco Choppers

Rosen Phaco Chopper

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM6194 45° angled, wedge shaped inferior edge, 1.3mm blunt tip.




3 in 1 Phaco Chopper

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM6153 Nucleus Chopper, Retractor and Nucleus Cracker in one 1.9mm tip



Kim Ball Tip Phaco Chopper

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM6154 Angled shaft, 2mm tip length
0.5mm ball tip, no cutting edge



Minardi Quick Chopper

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM6199 1.25mm Sharp tip



Premium Hooks and Manipulators

Kuglen Iris Hook and Lens Manipulator

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM5105 Angled Shaft




Kuglen Push/Pull Hook

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM5106 Angled Shaft



Sinsky Hook

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM5114 0.15mm diameter, blunt tip,
angled shaft



Bechert Nucleus Rotator

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM5153 Horizontal blunt forked tip

PM5154 Vertical blunt forked tip



Lester Lens Manipulator

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM5163 Angled, for lens with or without holes



Premium Hooks and Manipulators

Drysdale Nucleus Manipulator

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM5172 Paddle shaped tip




Mushroom Manipulator

Premium design, Titanium handle

Tip made of ultra-hardness Stainless Steel

PM5218 1.2mm angled shaft, mushroom tip



Koch Spatula

HK5236

Gently Curved 0.55mm width tip



HK5237

Angled shaft and angled tip



Premium Speculums



Adjustable Closed Wire Speculum

Premium design, Ultra-thin shaft with super light weight

Integrated forming body, elegant and sturdy

PM7344 Adjustable 15mm blades, wire closed,
78mm total length

Premium Markers

Mendez Degree Gauge

Premium design

made of ultra-hardness Stainless Steel

MR2804 For accurate axial alignment of the cornea, 10mm I.D. 16mm O.D. Calibrated every 5 degrees from 0 to 180 degrees. Well-defined groove scale.



Axis Marker

MR2805 This two line marker is used in conjunction with the Mendez Degree Gauge to indicate the axis of astigmatism.



Toric Marker


Premium design

made of ultra-hardness Stainless Steel

MR2806 This marker was designed for Toric Lens implant procedures to mark the lens and incision axis in one pass.



Premium Vitreo-Retinal Instruments



Ask our Sales Team to Get more Information about Vitreo-Retinal Instruments

Tying Forceps

McPherson Straight Tying Forceps

4.5mm tying platforms

TY1111 Standard tips, 85mm length

TY1113 Standard tips, 115mm length



McPherson 45° Angled Tying Forceps


4.5mm tying platforms

TY1121 Standard tips, 12mm shaft, 85mm length

TY1122 Standard tips, 12mm shaft, 115mm length

TY1123 Ultrafine tips, 7mm shaft, 85mm length

TY1124 Ultrafine tips, 10mm shaft, 85mm length



Strong Shaft Tying Forceps

Special stronger shaft, double size of the TY1113 & TY1122, 4.5mm tying Platforms to enhance gripping

TY1114 Straight shaft, 115mm length

TY1125 12mm Angled shaft, 115mm length



Curved Tying Forceps

6.0mm tying platforms

TY1131 85mm length

TY1132 115mm length



Castroviejo Tying Forceps

5mm shaft with tying platforms

TY1141 Straight, 108mm length

TY1142 Angled 45°, 108mm length



Tying Forceps

Castroviejo Tying Forceps

6mm tying platforms, made of Stainless Steel

TY1143 Straight, 108mm length

TY1144 Angled 45°, 10mm shaft , 108mm length



Girard Tying Forceps

4.5mm tying platforms

TY1151 Straight, 105mm length

TY1152 Curved, 105mm length



Toothed Forceps

Castroviejo Suturing Forceps

1 × 2 Forwarded Teeth, Straight Shaft with 5.0mm tying platforms

ST1421 0.12mm teeth, 108mm length

ST1422 0.3mm teeth, 108mm length

ST1423 0.5mm teeth, 108mm length



Short Shaft Toothed Forceps

Special short shaft gives more room when inserting into the wound

TS1322 0.12mm teeth, delicate tip, 108mm length



Toothed Forceps

Bonn Toothed Forceps

1 × 2 Standard Teeth, Straight Shaft with 5.0mm tying platforms

- FX1823 0.1mm teeth, 108mm length
- FX1821 0.12mm teeth, 108mm length
- FX1811 0.3mm teeth, 108mm length
- FX1822 0.4mm teeth, 108mm length
- FX1812 0.5mm teeth, 108mm length



Bishop-Harmon Tissue Forceps

Straight shaft without tying platform

- TS1311 Bonn Forceps 1 x 2 teeth, .12mm teeth, 85mm length
- TS1331 Straight, 1 x 2 teeth, 0.4mm teeth, 85mm length



Tissue Forceps

- TS1321 Straight, 0.3mm teeth, 100mm length



Tissue Forceps

Made of Stainless Steel

- TS1341 Straight 1 x 2 teeth, .1mm teeth, 108mm length
- TS1342 Straight 1 x 2 teeth, .12mm teeth, 108mm length
- TS1343 Angled 1 x 2 teeth, .12mm teeth, 108mm length



Straight Toothed Forceps

- ST1411 .12mm teeth, 85mm length
- ST1412 .10mm teeth, 85mm length
- ST1414 .12mm teeth, 105mm length
- ST1415 .10mm teeth, 105mm length



Toothed Forceps

Straight Toothed Forceps

ST1431 0.12mm ,1 x 2 teeth, 5.0mm tying platforms, 105mm long



Angled Kelman-McPherson Suturing Forceps

ST1441 .12mm 1 x 2 teeth, 5mm tying platforms, 10mm shafts, 85mm length



Curved Toothed Forceps

ST1451 .12mm teeth, 1 x 2 teeth, 5mm tying platforms



Fixation Toothed Forceps

FX1871 0.12mm 2x3 teeth, 105mm length

FX1872 0.3mm 2x3 teeth, 105mm length



Colibri Forceps

1 x 2 teeth

CN1211 0.12mm teeth with tying platforms, 85mm length

CN1212 0.12mm teeth with tying platforms, 115mm length

CN1213 45 degrees forwarded 0.12mm teeth with tying platforms, 85mm length



Toothed Forceps

Harms-Colibri Corneal Forceps

CN1221 1 x 2 teeth, .12mm teeth with tying platforms, 75mm length


CN1222 1 x 2 teeth, .30mm teeth with tying platforms, 75mm length

CN1223 1 x 2 teeth, .50mm teeth with tying platforms, 75mm length



Colibri Corneal Forceps

CN1231 1 x 2 teeth, .12mm teeth with tying platforms, 105mm length



Pierse-Colibri Corneal Forceps

CN1241 0.25mm notched, with tying platforms, 75mm length



St. Martin Corneal Forceps

CN1251 .12mm teeth, 5mm tying platforms, 85mm long

CN1252 .15mm teeth, 5mm tying platforms, 85mm long

CN1253 .3mm teeth, 5mm tying platforms, 85mm long



Polack Corneal Suturing Forceps

CN1261 1mm spread, delicate 1 x 2 .12mm teeth, 115mm length



Polack Corneal Suturing Forceps

CN1263 0.12mm, 1x2 teeth, Curved shafts without platforms, flat handle

Overall length 75mm



Toothed Forceps

Pierse Straight Fixation Forceps

with tying platforms

FX1831 Standard, 85mm length

FX1841 Standard, 115mm length



Pierse Angled Fixation Forceps

with tying platforms

FX1851 Standard, 85mm length



Curved Fixation Forceps with lock

FX1861 Curved shaft without platform, 0.5mm 1×2 teeth, overall length 118mm.



Moody Fixation Forceps, Stainless Steel

FX1862 Curved shaft with platform, 0.25mm 1×2 teeth with lock on left (first photo)

FX1863 Curved shaft with platform, 0.25mm 1×2 teeth with lock on right (second photo)

FX1864 Straight shaft with platform, 0.25mm 1×2 teeth with lock



Fixation Forceps, Stainless Steel

FX1865 4×4 teeth with lock on left (as photo shows)

FX1866 4×4 teeth with lock on right



T Fixation Forceps

FX1867 Angled 90° shafts with 10.0mm scaled jaws, 1mm teeth on one jaw which fits the holes on the other jaw, Overall length 110mm



Toothed Forceps


Fixation Forceps

FX1868 24mm "7" shaped shaft, with lock, Overall length 100mm



Adson Serrated Forceps

FX1879 0.8mm width serrated tip, Overall length 120mm



Adson Toothed Forceps

FX1880 0.6mm width 1×2 teeth tip, Overall length 120mm



Chalazion Forceps

Chalazion Forceps

CL1511 Round 8mm, OD 11mm, ID 8mm, 98mm length

CL1512 Round 15mm, OD 18mm, ID 15mm, 98mm length

CL1513 Round 12mm, OD 15mm, ID 12mm, 98mm length



Chalazion Forceps

CL1521 Oval 12 x 14mm, 98mm length

CL1522 Oval 9 x 17mm, 98mm length



Desmarres Chalazion Forceps

CL1531 inside dia.12x23mm



Meibomian Gland Forceps

Meibomian Gland Compressor

Broad, smoothly polished tips; designed to gently express meibum from inflamed or occluded meibomian glands, Made of Stainless Steel

MB1541 Round tip, 3.5mm Diameter, 105mm length. ↓



MB1542 Rectangle 3×3.5mm tip, 105mm length. ↓



MB1543 Rectangle 3×6mm tip, 105mm length. ↓



MB1544 Rectangle 5×6mm tip, 105mm length. ↓



MB1545 Rectangle 3.5×10mm tip, 105mm length. ↓



Meibum Expressor

The polished knurled rollers are designed to gently express meibum from inflamed or occluded meibomian glands

MB1546 Round knurled rollers 3.5mm Diameter, 105mm length.



Lens Forceps

Artisan Implantation Forceps

The implantation forceps is used for the implantation of the Toric ARTISAN IOL and fixation of the lens optic during the enclavation procedure.

LN1671 Implantation forceps refractive, long



Kelman-McPherson Forceps

LN1611 7.5mm smooth jaws, 85mm length

LN1612 7.5mm smooth jaws, 115mm length




Kelman-McPherson Lens Forceps

LN1621 Lens holding Forceps. Smooth jaws for holding soft IOLs, 102 mm long



Lens Implantation Forceps

LN1641 Smooth surface, round edges gently grip lens, 115mm length



Lens Loading Forceps

LN1652 for Alcon Monarch II



Lens Loading Forceps

LN1661 15mm bend, round smooth jaws for holding soft IOLs, 102mm


LN1662 12mm bend, round smooth jaws for holding soft IOLs, 102mm



Conjunctiva Forceps


Fechtner Conjunctiva Forceps

CJ1711 Very delicate ring tips with tying platforms, 112mm length




Moorfields Conjunctiva Forceps

CJ1721 1.0mm serrated, 100mm long



Dressing/Moorfields Forceps


CJ1731 fine 0.5mm serrated, 100mm long



Dressing/Moorfields Forceps Curved

CJ1722 Curved fine 1.0mm Serrated 100mm long

CJ1732 Curved fine 0.5mm Serrated 100mm long



Prechoppers

Akahoshi Combo Prechopper, Straight

PR1911 For nucleus grades 1-2. The sharp edge is used to insert into the nucleus, followed by the blunt edge to complete the division. 120mm length



Akahoshi Phaco Prechopper

PR1922 6mm long delicate jaws with sharp tips; for cracking, dividing and rotating the nucleus without sculpting
Made of Stainless Steel



Prechoppers

Nucleus Cracker

PR1923 Cross Action, Reverse control



Akahoshi Prechopper Cross-action Angled Shafts

PR1921 12mm shafts, maximum tip opening 3mm, 120mm length




Akahoshi Paddle Prechopper

PR1931 Straight shafts for 2.0mm incision, cross-action, maximum tip opening 2.5mm, 120mm long



Akahoshi Prechopper Cross-action Curved Shafts

PR1941 12mm shafts, maximum tip opening 3mm, 120mm length



Lens Cutters

IOL Cutting Scissors

One blade with serrations to securely grasp the IOL

LC3711 Straight blade.

LC3712 Angled blade.



IOL Cutting Scissors

LC3713 Designed to grasp and remove all types of foldable intraocular lens through a small incision, squeeze titanium Handle, overall length 120mm



Capsulorhexis Forceps

Corydon Style Capsulorhexis Forceps

CH1011 45° angled tips, curved 12 mm shafts, with guide pin and iris stop, 105 mm length



O'Gawa-Utrata Capsulorhexis Forceps, fine tips, cystotome teeth

CH1012 45° angled tips, straight 12mm shafts, with guide pin and iris stop, 105mm length



Utrata Capsulorhexis Forceps, cystotome teeth

CH1021 Shorter 10mm angled curved shaft, 105mm length



Capsulorhexis Forceps

CH1013 thinner long shaft for better vision, 12mm curved, 100mm long



Capsulorhexis Forceps

Made of Stainless Steel, 110mm total length

CH1014 12mm curved Shaft

CH1015 12mm straight Shaft



Capsulorhexis Forceps

Nevyas Capsulorhexis Forceps

108 mm, with iris stop

CH1031 45° angled tips, with 12 mm straight shafts

CH1032 45° angled tips, with 12 mm curved shafts



Nevyas Capsulorhexis Forceps

108mm total length with iris stop, Made of Stainless Steel

CH1033 45° angled tips, with 12 mm straight shafts

CH1034 45° angled tips, with 12 mm curved shafts



Utrata Capsulorhexis Forceps

fine sharp tips, cystotome teeth, 85mm length

CH1041 12mm angled and straight shaft

CH1042 12mm angled and curved shaft



Capsulorhexis Forceps, fine sharp tips, Stainless Steel

very delicate handle with thin long shaft

CH1043 14mm angled and straight shaft, 85mm length

CH1044 14mm angled and curved shaft, 85mm length



Capsulorhexis Forceps

Lehner-Utrata Capsulorhexis Forceps

very fine sharp smaller cystotome teeth, teeth smaller than CH1011

CH1051 12mm angled curved shaft, 105mm length



Inamura Capsulorhexis Forceps

1.8mm width at pivot box, Cross action tips, Round handle, 120mm

CH1061 Tip to pivot point 11mm

CH1067 Tip to pivot point 8.5mm



Inamura Capsulorhexis Forceps

CH1062 Tip to pivot point 8.5mm, front part 7mm shaft, cross action tips, 1.6mm width at pivot box, Flat handle, 120mm



Inamura Capsulorhexis Forceps

CH1063 Tip to pivot point 10.5mm, front part 7.5mm shaft, cross action tips, 1.5mm width at pivot box, Round handle, 120mm



Capsulorhexis Forceps

Inamura Capsulorhexis Forceps

CH1065 Tip to pivot point 12mm, there are 6 markings at the tip, 1mm per scale cross action tips, 1.5mm width at pivot box, Round handle, 120mm




Inamura Capsulorhexis Forceps

CH1066 Tip to pivot point 9mm, there are 5 markings at the tip, 1mm per scale cross action tips, 1.5mm width at pivot box, Round handle, 120mm



Giannetti MICS Capsulorhexis Forceps

CH1071 Ultra delicate shanks and tips with controlled opening, 10mm angled curved shaft, 105mm



Capsulorhexis Forceps for 1.8mm

Special thin long shaft for 1.8mm incision, made of ultra-hardness Stainless Steel

CH1072 14mm angled straight shaft, 110mm length


CH1073 14mm angled curved shaft, 110mm length



Miscellaneous Forceps

Cilia Forceps

CL1001/11412 72mm length



Cilia Forceps *Stainless Steel*

CL1009 90mm length



Towel Clamps

TW1002 72mm length



Serrefine

SF1003 10mm, Straight, 36mm long



Scleral Plug Forceps

PF1004 Securely hold the scleral plug during insertion and removal, 105mm long



Scleral Plug Forceps

PF1104 Double "C" shape rings to hold the scleral plug during insertion and removal, 118mm long



Jeweler Forceps

JL1006 .05 x .01mm Straight tip

JL1007 .05 x .01mm Curved tip

JL1008 .05 x .01mm Angled tip



Miscellaneous Forceps

Green Fixation Forceps

FX1010 Extured grip atraumatic fine toothed jaws to allow ease of access to the globe.



Grooved Inserter Forceps

GF1090 For insertion of extended duration and collagen implants



Corneal Transplant Forceps

CN1264 T-shaped forceps for avulsion of posterior elastic layer during corneal transplant.



Corneal Transplant Forceps

CN1265 Ring forceps for corneal transplant.



Iris Scissors

Straight Iris Scissors

IR3111 11mm blades, sharp tips, 4mm cutting length, 90mm length




Curved Iris Scissors

IR3112 11mm blades, sharp tips, 90mm length



Straight Iris Scissors

IR3121 11mm blades, sharp tips, 105mm length



Curved Iris Scissors

IR3122 11mm blades, sharp tips, 105mm length



Corneal Scissors

Castroviejo Corneal Section Scissors - Left Hand

CN3212 Blunt tips, lower blade .5mm longer than upper, 9mm from pivot to tip, 105mm length

CN3222 Blunt tips, lower blade .5mm longer than upper, 13mm from pivot to tip, 105mm length



Castroviejo Corneal Section Scissors - Right Hand

CN3211 Blunt tips, lower blade .5mm longer than upper, 9mm from pivot to tip, 105mm length

CN3221 Blunt tips, lower blade .5mm longer than upper, 13mm from pivot to tip, 105mm length



Castroviejo Curved Corneal Scissors

CN3231 Blunt tips, 16mm from pivot to tip, 115mm length



Castroviejo Curved Corneal Scissors

CN3240 Blunt tips, small blades, 17mm from pivot to tip, 105mm length



Castroviejo Curved Corneal Scissors

CN3241 Blunt tips, 11mm from pivot to tip, 105mm length



Castroviejo Straight Corneal Scissors

CN3242 Blunt tips, 11mm from pivot to tip, 105mm length



Corneal Scissors

Castroviejo Keratoplasty Scissors

CN3251 Blunt angled tips, longer lower blade, 16mm from pivot to tip
105mm length



Castroviejo Curved Corneal Scissors

CN3261 Curved blunt tips, 10.5mm from pivot to tip, 105mm length



Castroviejo Curved Corneal Scissors, blunt tips

CN3271 16mm from pivot to tip, 115mm length


CN3272 11mm from pivot to tip, 112mm length



Castroviejo Curved Corneal Scissors

Made of Stainless Steel

CN3232 Blunt tips, 16mm from pivot to tip, 115mm length



Corneal Scissors

Castroviejo Curved Corneal Scissors *Stainless Steel*

12mm from pivot to tip, 120mm length

CN3281 blunt tips

CN3282 sharp tips



Castroviejo Curved Corneal Scissors *Stainless Steel*

14mm from pivot to tip, 120mm length

CN3283 blunt tips

CN3284 sharp tips



Castroviejo Curved Corneal Scissors *Stainless Steel*

16mm from pivot to tip, 120mm length

CN3285 blunt tips



Tenotomy & Utility Scissors

Westcott Curved Tenotomy Scissors

TN3311 Blunt tips, 19mm from pivot to tip, 110mm length

TN3321 Blunt tips, 21mm from pivot to tip, 115mm length



Westcott Utility Scissors

TN3331 Sharp tips, 25mm from pivot to tip, curved blades, 125mm length



Utility Scissors

UT3341 Straight blades, pointed tips, 32mm from pivot to tip, 110mm length

UT3342 Straight blades, blunt tips, 32mm from pivot to tip, 110mm length

UT3351 Curved blades, pointed tips, 32mm from pivot to tip, 110mm length

UT3352 Curved blades, blunt tips, 32mm from pivot to tip, 110mm length



Enucleation Scissors Stainless Steel

EN3354 Curved blunt tips, 52mm from pivot to tip, 130mm length



Tenotomy & Utility Scissors

Utility Scissors *Stainless Steel*

Standard stainless steel scissors, good for limited budget.

UT3361 Straight blades, pointed tips, 28mm from pivot to tip, 100mm length

UT3362 Curved blades, pointed tips, 28mm from pivot to tip, 100mm length



Utility Scissors *Stainless Steel with coated ring*

Coated with special layer, the blades extraordinary sharp and last much longer.

UT3365 Straight blades, pointed tips, 28mm from pivot to tip, 95mm length

UT3366 Straight blades, blunt tips, 28mm from pivot to tip, 95mm length

UT3367 Curved blades, pointed tips, 28mm from pivot to tip, 95mm length

UT3368 Curved blades, blunt tips, 28mm from pivot to tip, 95mm length



Vannas Scissors

Vannas Straight Scissors

VN3411 Sharp pointed tips, 7.5mm blades, 9.5mm from pivot to tip, 82mm



Vannas Curved Scissors

VN3412 Sharp pointed tips, 7.5mm blades, 9.5mm from pivot to tip 82mm



Vannas Angled Scissors

VN3421 Sharp tips, 7.0mm blades, 9.0mm from pivot to tip 82mm



Vannas Scissors

Vannas Straight Scissors *Stainless Steel*

VN3413 Sharp pointed tips, 8.5mm blades,
85mm total length



Vannas Curved Scissors *Stainless Steel*

VN3414 Sharp pointed tips, 8.5mm blades,
85mm total length



Vannas Angled Scissors *Stainless Steel*

VN3422 Sharp tips, 8.0mm blades,
85mm total length



Vannas Straight Capsulotomy Scissors

VN3431 Sharp pointed tips, 10.5mm blades, 108mm length

VN3432 Sharp pointed tips, 12.5mm blades, 108mm length



Vannas Straight Capsulotomy Scissors *Stainless Steel*

Sharp pointed tips

VN3433 12mm from pivot to tip, 110mm length

VN3434 16mm from pivot to tip, 110mm length



Vannas Curved Capsulotomy Scissors

VN3441 Sharp pointed tips, 10.5mm blades, 95mm length

VN3442 Sharp pointed tips, 12.5mm blades, 95mm length



Vannas Scissors

Vannas Curved Capsulotomy Scissors

Sharp pointed tips, made of Stainless Steel

VN3443 12mm from pivot to tip, 110mm length



Vannas Angled Capsulotomy Scissors

VN3451 Sharp pointed tips, 10.5mm blades, 95mm length

VN3452 Sharp pointed tips, 12.5mm blades, 95mm length




Vannas Angled Capsulotomy Scissors Stainless Steel

VN3453 Sharp pointed tips, 16mm from pivot to tip, 110mm length



Vannas Scissors

VN3461 Extra thin 7mm blades, sharp pointed tips, 9mm from pivot to tip 105mm lengths




Vannas Scissors

Extra thin blades, sharp point, 15mm from pivot to tip, 108mm long

VN3471 Straight blade

VN3472 Curved blade



Suture & Stitch Scissors

Westcott Stitch Scissors

SS3511 Very sharp pointed tips, 19mm from pivot to tip, 115mm length



Mcperson-Westcott Stitch Scissors

SS3521 Very sharp pointed tips, 10.5mm from pivot to tip, 105mm length



Mcperson-Westcott Stitch Scissors

SS3531 Sharp pointed tips, 13mm from pivot to tip, 115mm length
SS3541 Sharp pointed tips, 13mm from pivot to tip, 140mm length



Conjunctiva Scissors

Conjunctiva Scissors

CJ3611 Gently curved blades, tip to pivot length 12mm, blunt tips, 115mm long



Westcott Conjunctival Scissors

CJ3621 Lightly curved blunt tips, 16mm from mid-screw to tip, 130mm long



Osher Universal Scissors

CJ3631 Blunt tips, thin curved blades, 21mm from pivot to tip, 115mm length




Needle Holders

Troutman Curved Needle Holder

NH4111 Delicate 10mm jaws, blunt, without lock, 115mm length


NH4112 Delicate 10mm jaws, blunt, with lock, 115mm length



Barraquer Curved Needle Holder

NH4121 Delicate 9mm jaws, blunt, without lock, 115mm length

NH4122 Delicate 9mm jaws, blunt, with lock, 115mm length



Fine Curved Needle Holder

NH4131 Fine 9mm jaws, blunt, without lock, 115mm length

NH4132 Fine 9mm jaws, blunt, with lock, 115mm length



Barraquer Curved Needle Holder

NH4141 Standard 9mm jaws, blunt, without lock, 115mm length


NH4142 Standard 9mm jaws, blunt, with lock, 115mm length



Barraquer Straight Needle Holder

NH4151 Standard 9mm jaws, blunt, without lock, 115mm length


NH4152 Standard 9mm jaws, blunt, with lock, 115mm length



Curved Needle Holder

NH4161 Standard 11mm jaws, blunt, without lock, 115mm length

NH4162 Standard 11mm jaws, blunt, with lock, 115mm length



Needle Holders

Straight Needle Holder

NH4171 Standard 6mm jaws, without lock, 110mm length

NH4172 Standard 6mm jaws, with lock, 110mm length



Barraquer Straight Needle Holder

NH4181 Standard 9mm jaws, without lock, 115mm length



Bechert-Sinskey Curved Needle Holder

NH4182 Standard 9mm jaws, without lock, 115mm length



Bechert-Sinskey Needle Holder, Very Fine tips

Very Fine 9mm gently curved jaws, sharp tips, 115mm long

NH4183 Without Lock

NH4184 With Lock



Curved Needle Holder

NH4191 Very Fine 6mm jaws, sharp tips, without lock, 110mm length

NH4192 Very Fine 6mm jaws, sharp tips, with lock, 110mm length




Micro Needle Holders, curved

Very Fine 6mm moderately curved jaws, without lock, 128mm long

NH4193 Without Lock

NH4194 With Lock



Needle Holders

Barraquer Curved Needle Holders

Delicate 9mm curved jaws, 125mm long

NH4101 Without Lock

NH4102 With Lock



Castroviejo Curved Needle Holder

Gently curved jaws, Strong blunt tips, 140mm long

NH4211 Standard 11mm jaws, blunt, without lock, 140mm length

NH4212 Standard 11mm jaws, blunt, with lock, 140mm length



Castroviejo Straight Needle Holder (Closed tip width 0.8mm)

NH4221 Standard 9mm jaws, without lock, 125mm length

NH4222 Standard 9mm jaws, with lock, 125mm length



Castroviejo Curved Needle Holder (Closed tip width 0.8mm)

NH4231 Standard 9mm jaws, without lock, 125mm length

NH4232 Standard 9mm jaws, with lock, 125mm length



Azar Curved Needle Holder (Closed tip width 0.7mm)

NH4241 Standard 9mm jaws, without lock, 115mm length

NH4242 Standard 9mm jaws, with lock, 115mm length



Needle Holders

Castroviejo Needle Holder

Gently curved jaws, Fine blunt tips, 140mm long

NH4251 Without Lock

NH4252 With Lock



Castroviejo Straight Needle Holder (Closed tip width 1.0mm)

NH4261 Standard 11mm jaws, blunt, without lock, 140mm length

NH4262 Standard 11mm jaws, blunt, with lock, 140mm length



Hooks Spatulas & Manipulators

Kuglen Iris Hook and Lens Manipulator

HK5101 Angled

HK5102 Straight



Kuglen Push/Pull Hook

HK5103 Angled

HK5104 Straight



Sinskey Hook

HK5111 0.15mm diameter, blunt tip, angled shaft

HK5112 0.15mm diameter, blunt tip, straight shaft



Hooks Spatulas & Manipulators

Reverse Sinskey Hook

HK5113 .15mm blunt tip, angled shaft



Fenzel Insertion Hook

HK5121 Angled

HK5122 Straight



Jaffe-Knolle Iris Hook

HK5131 Angled, curved blunt hook

HK5132 Straight, curved blunt hook



Lewicky Lens Manipulating Hook

HK5141 90° angled tip on curved shaft, .15mm blunt tip



Bechert Nucleus Rotator

HK5151 Horizontal blunt forked tip

HK5152 Vertical blunt forked tip



Lester Lens Manipulator

HK5161 Angled, for lens with or without holes

HK5162 Straight, for lens with or without holes



Drysdale Nucleus Manipulator

HK5171 Paddle shaped tip



Hooks Spatulas & Manipulators

Button Manipulator

HK5181 Angled 45° shaft, roughened posterior surface, designed for iris retraction and capsule polishing.



Maltzman-Fenzel Lens Manipulator

HK5191 Angled, .15mm V-shaped tip

HK5192 Straight, .15mm V-shaped tip



Connor Wand, curved tip

Designed to engage the nucleus equator to stabilize it and provide counter pressure while the prechopper forceps is being inserted into the nucleus core.

HK5201 Angled shaft

HK5202 Straight shaft



Nucleus Manipulator

HK5211 Angled shaft with curved tip

Designed to safely spin, chop, and Back-Crack a nucleus using less phaco time. Perfect for dialing intraocular lenses in and out of the bag



90° Nucleus Manipulator

HK5212

Designed to safely spin, chop, and back crack a nucleus using less phaco time. Used for sweeping nucleus fragment toward a phaco tip in the central safety zone for safe emulsification, with less stress to the zonules. Perfect for dialing intraocular lenses in and out of the bag.



Hooks Spatulas & Manipulators

Small Pupil Snapper Hook

HK5213 Smooth polished notch superiorly to push back iris with sharp inner edge for nucleus splitting.



Fukasaku Small Pupil Snapper Hook

HK5216 Smooth polished notch superiorly to push back iris with sharp inner edge for nucleus splitting.



Escaf Nucleus Manipulator

HK5214 Flat oval tip, angled shaft



Oliver Nucleus Manipulator

HK5215 Angled, with straight, olive-shaped tip



ICL Manipulator

HK5221 Manipulates inferior and superior haptic of lens for positioning under the iris. Posterior surface is micro sandblasted.



Mushroom Manipulator

HK5217 1.2mm angled shaft, mushroom tip



Hooks Spatulas & Manipulators

Barraquer Cyclodialysis Spatula

- HK5231 Semi-sharp .5mm tips 0.1mm thick, Curved
- HK5232 Semi-sharp .7mm tips 0.1mm thick, Curved
- HK5234 Semi-sharp .5mm tips 0.1mm thick, Straight
- HK5235 Semi-sharp .7mm tips 0.1mm thick, Straight



Mackool Barraquer Spatula

- HK5233 spatula with a polished tip 0.5mm wide and 0.25mm thick



Barraquer Cyclodialysis Spatula

- HK5241 Semi-sharp .25mm vaulted tips. Double-ended 10mm tip to bend and 15mm tip to bend



Double end Spatula

- HK5242 Double end spatula, gently curved, 1.5mm width and 2mm width.



Minami "M" Hook

- CP6121 /14127 Chopper, Splitter, Dissector, Spatula, Rotator and Manipulator in one.



Strabismus Hooks

Helveston Muscle Hook

- HK5001 Delicate, blunt 8mm spatulated tip
- HK5002 Delicate, blunt 10mm spatulated tip
- HK5003 Delicate, blunt 12mm spatulated tip



Retinal Detachment Hook

- HK5004 Hook with a ring at the tip



Helveston Finder Hook

- HK5005 Lightly angled cone shaped tip



Stevens Curved Tenotomy Hook

- HK5006 6mm curved hook, 120mm



Graefe Strabismus Hook

- HK5007 Size 8mm*1.5mm, 135mm
- HK5008 Size 9mm*1.5mm, 135mm
- HK5009 Size 10mm*1.5mm, 135mm
- HK5010 Size 10mm*2.0mm, 135mm



Retractors

Helveston Retractors

- RT5016 7mm Retractor
- RT5017 9mm Retractor
- RT5018 11mm Retractor



Retractors

Desmarres Lid Retractor

- RT5019 8mm wide Retractor
- RT5020 11mm wide Retractor
- RT5021 13mm wide Retractor
- RT5022 15mm wide Retractor



Large Retractors

- RT5026 6mm×25mm wide, flat handle overall length 135mm
- RT5027 8mm×25mm wide, flat handle overall length 135mm



Schepens Orbital Retractor

- RT5023 14.0mm x 57.0mm blade, smooth with notched end, 4.5mm wide notch, 150mm long



Knapp (Rollet) Lacrimal Sac Retractor

- RT5024 /34049 4 x 8mm wide, sharp, 4-prong tip
- RT5025 /34087 4 x 8mm wide, blunt, 4-prong tip



Iris Retractors

- RT5031 Reusable by EO disinfection, made of PMMA, 0.2mm diameter, 1 set include 5 pcs, 13mm length

Iris Retractors Sterilization Caddy

- RT5032 Made of PTFE, designed to place iris hook, can be used for 5pcs iris retractors, Size Φ 30x21mm



Phaco Choppers



Phaco Chopper and Manipulator
CP6101 /14302 Double ended Drysdale Manipulator & Nagahara Phaco Chopper



Phaco Chopper and Manipulator
CP6102 /14303 Double ended Lieberman Micro Finger & Sharp Quick Minardi Chopper



Phaco Chopper and Manipulator
CP6103 /34128 Phaco Chopper with mushroom manipulator



Seibel Nucleus Chopper and Safety Quick Chopper
CP6104 Double ended, angled shafts. One end is olive-shaped tip with curved interior cutting edge. The other end is olive-shaped blade. Round handle.



Nagahara Phaco Chopper
CP6111 /34024 Curved shaft, 1.1mm tip
CP6112 /34025 Curved shaft, 1.6mm tip



Phaco Chopper
CP6161 /14122 Angled shaft
CP6162 /14121 Straight shaft



Claw Phaco Chopper
CP6131 /14130 Special triple inner blades construction allows chopping in three directions.



Phaco Chopper
CP6141 /14141 Special inner blade offset 45°
Unique elongated distal segment



Haefliger Phaco Cleaver
CP6171 /14139 Chopper and manipulator in one



Seibel Nucleus Chopper
CP6191 Round Polished tip



Rosen Phaco Splitter
CP6193 /34022 45° angled, wedge shaped inferior edge, 1.3mm blunt tip



Minardi Quick Chopper
CP6195 Sharp quick chopper



Straight



Angled Straight




Angled Curved

Phaco Choppers

Phaco Chopper & Manipulator

CP6101 Paddle Manipulator/ Phaco Chopper



Phaco Chopper & Manipulator

CP6102 Lieberman Micro Finger/Quick Chop



Phaco Chopper & Manipulator

CP6103 Phaco Chopper with mushroom manipulator




Double end Seibel Chopper *Stainless Steel*

CP6104 Double-ended, angled shafts. Overall length: 138mm.

One end is angled 15mm from bend to tip with curved interior cutting edge and 0.5mm olive-shaped tip.

The other end is angled 10mm from bend to tip with 0.5mm olive-shaped blade and 1.0mm bend.




Phaco Chopper & Manipulator

CP6105 2mm Nagahara Chopper with Ball tip Manipulator



Nagahara Phaco Chopper

CP6106 1.2mm Nagahara Chopper and 1.9mm Nagahara Chopper



Phaco Choppers


Sinskey Hook and Drysdale Manipulator

CP6107 Sinskey Hook and Drysdale Manipulator



Cataract Coach Chopper

CP6108 This universal cataract instrument is designed to allow multiple functions in one simple and efficient instrument. The chopping tip is designed to provide excellent control of the nucleus with vertical, horizontal or tilt-chop, as well as other techniques. The gently curved and smooth tip, designed to protect the posterior capsule, features grooves intended to engage the optic or haptics for IOL manipulation and for iris push-pull in small pupil cases. A traditional spatula is designed for intraocular manipulation of iris, lens fragments and vitreous. Overall length 137mm.



Nagahara and Ball tip Phaco Chopper

CP6109 2.0mm Nagahara Chopper and 1.5mm Ball tip



1.5 mm Ball Tip

2.0 mm Nagahara Tip



Kim Ball Tip Phaco Chopper


CP6151 Angled 2mm shaft, 2mm diameter ball tip, no cutting edge



Phaco Choppers

3 in 1 Phaco Chopper

CP6152 Nucleus Chopper, Retractor and Nucleus Cracker in one 1.9mm tip



Nagahara Phaco Chopper

CP6111 Angled Curved shaft, 1.1mm tip

CP6112 Angled Curved shaft, 1.6mm tip



Phaco Chopper

CP6161 Angled shaft

CP6162 Straight shaft

Nichamin Phaco Chopper

Rounded and highly polished distal tip.

CP6181 Straight, ambidextrous nucleus chopper, tip 1.1 mm long

CP6182 Curved, ambidextrous nucleus chopper, tip 1.1 mm long

CP6187 Straight, ambidextrous nucleus chopper, tip 1.8 mm long

CP6188 Curved, ambidextrous nucleus chopper, tip 1.8 mm long

Phaco Side Cutter Chopper

Phaco side cutter chopper, Round and highly polished distal tip.

CP6183 Angled-Straight Shaft, tip 1.25mm long, right hand

CP6184 Angled-Straight Shaft, tip 1.25mm long, left hand

CP6185 Angled-Curved Shaft, tip 1.25mm long, right hand

CP6186 Angled-Curved Shaft, tip 1.25mm long, left hand

CP6196 Angled-Straight Shaft, tip 1.6mm long, right hand

CP6197 Angled-Straight Shaft, tip 1.6mm long, left hand

CP6189 Angled-Curved Shaft, tip 1.8mm long, right hand

CP6180 Angled-Curved Shaft, tip 1.8mm long, left hand

Phaco Choppers

Claw Phaco Chopper

CP6131 Special triple inner blades construction allows chopping in three directions.



Phaco Chopper

CP6141 Special inner blade is offset 45°. Unique elongated distal segment



Haefliger Phaco Cleaver

CP6171 Chopper and manipulator in one



Seibel Nucleus Chopper

CP6191 Chopper with ball-shaped tip for efficient capsule protection



Rosen Phaco Chopper

CP6193 45° angled, wedge shaped inferior edge, 1.3mm blunt tip.




Minardi Quick Chopper

CP6195 Sharp quick Chopper



Green's Phaco Chopper



CP6198 Angled.



Speculums

Lieberman Aspirating Speculum-Stainless Steel

SP7002 For removal of excess fluid from the ocular surface. Reversible speculum can rest nasally or temporally. Adjustable 14.5mm blades, 65mm length




Wire Speculum

SP7111 13.5mm blades, 52mm length



Kratz-Barraquer Wire Speculum

SP7121 13.5mm open blades, 52mm length



Nagahara Wire Speculum (Nasal placement)

SP7131 Designed for temporal approach cataract surgery
13.5mm blades, 52mm length



Grandon Wire Speculum

SP7141 Gently curved back tips retain redundant lid tissue, heavy gauge coiled fulcrum provides strong spring tension. 13.5mm blades, 52mm length



Slade LASIK Eye Speculum

SP7151 13.5mm blades, 52mm length



Pediatric Wire Speculum *Stainless Steel*

SP7161 9mm blades, 29mm length

SP7162 7mm blades, 29mm length

SP7163 4mm blades, 29mm length



Pediatric Speculum *Stainless Steel*

SP7164 4mm blades, 35mm length


SP7165 11mm blades, 40mm length

Speculums




Wire Speculum

SP7171 7mm blades, 47mm length



Barraquer Solid Blade Speculum

SP7211 14.5mm Straight blades, 45mm length



Barraquer Solid Blade Speculum

SP7215 11mm Gently Curved blades



Kownacki LASIK Speculum

SP7422 14.5mm Gently Curved blades

No interfere with trying to pull down the lid to tighten up the conjunctiva during Femto LASIK.
Perfect room for the suction device.



Barraquer Solid Blade Speculum *Stainless Steel*

SP7212 4mm blades, 30mm length Pediatric

SP7213 10mm blades, 45mm length Pediatric

SP7214 15mm blades, 50mm length



Cook Adjustable Blade Speculum, *Temperal Side*

Stainless Steel with locking screw

SP7012 Adjustable 8mm blades, 40mm length

SP7013 Adjustable 11mm blades, 43mm length

SP7014 Adjustable 14.5mm blades, 45mm length

Speculums



Cook Adjustable Blade Speculum, Nasal Side

Stainless Steel with locking screw

SP7015	Adjustable	8mm blades, 40mm length
SP7016	Adjustable	11mm blades, 43mm length
SP7017	Adjustable	14.5mm blades, 45mm length



Clark Open Wire Speculum

SP7011 14mm blades, 82mm length



Lieberman Open V-Wire Speculum

SP7301 Adjustable 14mm blades, temporal side, 85mm length

SP7303 Adjustable 14mm blades, nasal side, 85mm length



Slade LASIK Wire Speculum

SP7302 Designed for maximum exposure of globe and to accommodate suction ring. Adjustable 14mm blades, 85mm length



Lieberman-Kratz Open Wire Speculum

SP7310 14mm open blades, temporal side, 45° angled shaft 80mm length



Lieberman-Kratz Open Wire Speculum

SP7311 16.5mm open blades, temporal side, 45° angled shaft, 85mm length

SP7312 16.5mm open blades, nasal side, 45° angled shaft, 85mm length



Adjustable Speculum *Pediatric*

SP7313	5mm wire
SP7314	7mm wire
SP7315	9mm wire
SP7316	11mm wire

Speculums



Reversible Lid Speculum

SP7321 Can be placed nasally or temporally to accommodate the suction ring of a microkeratome. Adjustable 12.5mm blades, 85mm length



Shapira Eye Speculum

SP7331 Adjustable wire opening for phaco and micro lamellar keratoplasty. Adjustable 13mm blades, 85mm length




Adjustable Closed Wire Speculum

SP7341 Adjustable 15mm blades, wire closed, 85mm length



Adjustable Closed Wire Speculum

SP7342 Adjustable 14.5mm blades, wire closed, 85mm length



Adjustable Closed Wire Speculum

SP7343 Adjustable 14mm blades, wire closed, 85mm length



Lieberman Solid Speculum

SP7411 Adjustable 14.5mm blades, 80mm length



Steinert LASIK ACS Speculum, temporal fit

SP7421 Strong with ideal curvature for accommodating the microkeratome suction ring. Adjustable 14.5mm blades 80mm length



Genisi LASIK Speculum

SP7431 Strong with ideal curvature for accommodating the microkeratome suction ring. Adjustable 14.5mm blades, 80mm length

Speculums



Machat Adjustable Speculum

SP7501 Provides maximum exposure of the globe for superior LASIK procedure. Open blade design avoids interference from speculum with the microkeratome, 80mm

ICL® Instruments

ICL® is a registered brand of Staar Surgical.


ICL® Manipulator

IL7801 Manipulates inferior and superior haptic of lens for positioning under the iris.



ICL® Manipulator

IL7802 Double-ended micro ring and hemi-ball tip, for atraumatic traction on the surface of an ICL®.



ICL® Olive Manipulator

IL7803 Designed to tuck haptics of the ICL® under the iris plane. Olive atraumatic tip engages the haptics allowing for easy placement of the ICL®.



ICL® Cartridge Loading Forceps

LN1662 12mm long, smoothly polished jaws, for careful placement of the ICL® into the injector cartridge.



ICL® Lens Cartridge Loading Forceps

IL7805 Straight 20G, jaws feature a series of concave and convex surfaces. For positioning the ICL in the cartridge without damaging the optic.



LASIK LASEK & ReLEx SMILE



Smile Femtosecond Forceps

LF7731 Curved textured tip, Designed for ReLEx® Smile® procedure




ReLEx Dissector & Lifter

LA7732 Spatula with semi-sharp spoon-shaped tip for dissection of lenticule.



ReLEx Dissector & Lifter

LA7733 Dissector & Lifter



Dissector & Lifter

LA7704 Dissector & Lifter

LASIK LASEK & ReLEx SMILE



Lasik Spatula
LA7703 Lasik Spatula



Corneal Epithelium Spatula
EP5261 Epithelium Spatula



Corneal Epithelium Spatula
EP5262 Semi-sharp end of the instrument designed to remove epithelium.




LASIK Flap Manipulator
LA7701 Olive shape tip to locate the flap edge and open the flap



Flap Elevator
LA7702 1.5mm wide, 12mm long, semi-sharp edges




LASIK LASEK & ReLEx SMILE



Femto Flap Lifter & Retreatment Spatula

Double ended instrument with a delicate pick and a curved spatula to open and manipulate the flap.



LA7705 7-shape + 135 degrees spatula



Femto Flap Lifter

Double ended instrument with a delicate pick and a curved spatula to open and manipulate the flap.



LA7735 7-shape + 90 degrees spatula



Femto Flap Lifter

Double ended instrument with a delicate pick and a curved spatula to open and manipulate the flap.


LA7734 Strong curved +135 degrees spatula



Femto Flap Lifter

Double ended instrument with a delicate pick and a curved spatula to open and manipulate the flap.

LA7736 Slightly curved +135 degrees spatula



Femto Flap Lifter

Double ended instrument with a delicate pick and a curved spatula to open and manipulate the flap.

LA7737 Slightly curved +90 degrees spatula

LASIK LASEK & ReLEx SMILE



Mendez Multipurpose LASIK Forceps

LA7706 Combined features of Spatula, Dissector and Forceps, 115mm long



Perone LASIK Flap Forceps

LA7707 2.0mm diameter rings tip, 108mm long




Manche LASIK Forceps

LA7708 Solid oval tip, 108mm long

Camellin Style LASEK Alcohol Well

AW2213	9.5mm diameter
AW2214	9.0mm diameter
AW2215	8.5mm diameter
AW2216	8.0mm diameter
AW2217	7.5mm diameter



Spatulated LASIK Cannula

LA7709 Flat spatulated tip, end opening 26G



Gimbel LASIK Fountain Cannula

LA7710 Angled 10mm bend tip, formed to the curvature of the cornea, 25G



Buratto LASIK Cannula

Flat smoothly burnished tip, easy to insert under the corneal flap. Two large side ports.

LA7711	25G
LA7712	23G



LASIK LASEK & ReLEx SMILE

LASIK Corneal Marker

MK2212 3 blades



Paufique Knife

PK2914 Corneal dissector, 6mm



Markers

Thornton Optical Zone Marker

With needle point in the center

MK2209 Optical Zone Marker, 5.0mm

MK2200 Optical Zone Marker, 5.5mm

MK2201 Optical Zone Marker, 6.0mm

MK2202 Optical Zone Marker, 6.5mm

MK2203 Optical Zone Marker, 7.0mm

MK2204 Optical Zone Marker, 7.5mm

MK2205 Optical Zone Marker, 8.0mm

MK2206 Optical Zone Marker, 8.5mm

MK2207 Optical Zone Marker, 9.0mm

MK2208 Optical Zone Marker, 9.5mm



Osher-Neumann 8-line Corneal Marker

MK2210 8-line blade, marks radial lines

4mm long, 4.5mm inner diameter of ring, 110mm long



Osher-Neumann 16-line Corneal Marker

MK2211 16-line blade, marks radial lines

4mm long, 4.5mm inner diameter of ring, 110mm long



Markers

Gravity Marker

Pre-Op Reference Marker

Designed for Limbal Relaxing Incisions and Toric IOL Implantation. Gravity weight system Assures excellent visualization and balance.

MK2221 Angled shaft

MK2222 Straight shafts



Bubble Toric Marker

MR2224 This marker was designed for Toric Lens implant procedures to mark the lens and incision axis in one pass. The instrument features a bubble level to aid in proper orientation of the marker. Total length 115mm.



Mendez Degree Gauge

MR2801 For accurate axial alignment of the cornea, 12mm I.D. 19mm O.D. Calibrated every 10 degrees from 0 to 180 degrees.

MR2802 For accurate axial alignment of the cornea, 12mm I.D. 16mm O.D. Calibrated every 10 degrees from 0 to 180 degrees.



Mendez Degree Gauge

MR2803 For accurate axial alignment of the cornea, 12mm I.D. 16mm O.D. Calibrated every 10 degrees from 0 to 180 degrees.



Graether Half-circle Axis Marker

MR2811 Can be used as a pre-op and intra-op marker for Toric. Inner tip-to-tip measures a 6 mm guide for capsulorhexis.



Markers


Axis Marker

MR2812 Bores Axis Marker This two line marker is used in conjunction with the Mendez Degree Gauge to indicate the axis of astigmatism.



LRI Marker

MR2821 Intra-Op LRI, Automatically creates marks at 40-60-80 degrees of Limbal Relaxing Incision.



LRI Marker/Toric IOL Marker

MR2831 Intra-Op LRI/Toric IOL Marker. Automatically creates marks for 6 & 8mm Chord length of Limbal Relaxing Incision.



3 Blade Marker for Clear

MR2807 This 3 Blade Marker is a precision-crafted specialty surgical instrument designed for ophthalmic procedures, particularly with Ziemer's CLEAR platform.



Irrigation & Aspiration

I/A Handpiece 21G, 0.3mm asp. Port

Titanium Handle with Stainless steel tip.

IA2301 I/A Handpiece 21G Straight Tips

IA2302 I/A Handpiece 21G Curved Tips

IA2303 I/A Handpiece 21G Angled Tips



I/A Handpiece for 2.2mm incision

Whole titanium design, use together with silicone sleeves.

IA2343 I/A Handpiece for 2.2mm incision Straight Tips

IA2344 I/A Handpiece for 2.2mm incision Curved Tips

IA2345 I/A Handpiece for 2.2mm incision 45 degree Angled Tips



I/A Handpiece for 2.8mm incision

Whole titanium design, use together with silicone sleeves.

IA2346 I/A Handpiece for 2.8mm incision Straight Tips

IA2347 I/A Handpiece for 2.8mm incision Curved Tips

IA2348 I/A Handpiece for 2.8mm incision 45 degree Angled Tips



I/A Handpiece Tips 21G, 0.3mm asp. Port

IA2304 I/A Tips Straight Tips Stainless Steel(for IA2301)

IA2305 I/A Tips Curved Tips Stainless Steel(for IA2302)

IA2306 I/A Tips Angled Tips Stainless Steel(for IA2303)

IA2307 **I/A Silicon Sleeves** for IA2301, IA2302 and IA2303



Infusion Handpiece

IA2308 7mm Handle, 90mm length, one end.

IA2309 7mm Handle, 95mm length, double end.



Irrigation & Aspiration

Backflush Handle

IA2310 overall length 158mm Passive Aspiration 20G

IA2336 overall length 158mm Passive Aspiration 23G



Backflush Handle

IA2337 overall length 158mm Active Aspiration 20G

IA2338 overall length 158mm Active Aspiration 23G



Bi-Manual Irrigation Handpieces

IA2311 Curved 21G, 2x0.4mm side ports.

IA2313 Curved 23G, 2x0.3mm side ports.

IA2327 Curved 21G, 1x0.5mm front port



Bi-Manual Aspiration Handpieces


IA2312 Curved Bi-Manual Aspiration Handpiece 21G, 1x0.3mm side-port

IA2314 Curved Bi-Manual Aspiration Handpiece 23G, 1x0.3mm side-port



Aspiration Handpiece with capsule polisher

IA2329 21G, 0.3mm side-port, with sandblasted capsule polisher



Irrigation & Aspiration

Bi-Manual Irrigation Handpieces

Improved smoother cannula tube with highly polished larger side ports.
Suitable for larger irrigation flow.

- IA2315 Curved Bi-Manual Irrigation Handpiece 20G, 2×0.3mm side-port
- IA2317 Curved Bi-Manual Irrigation Handpiece 23G, 2×0.3mm side-port



Bi-Manual Aspiration Handpieces

Improved smoother cannula tube with highly polished aspiration port.

- IA2316 Curved Bi-Manual Aspiration Handpiece 20G, 1×0.3mm port
- IA2318 Curved Bi-Manual Irrigation Handpiece 23G, 1×0.3mm port




Cannula I/A

Simcoe Double-Barreled I/A Unit

Thin wall irrigation and aspiration, aspiration through 0.3mm top port with irrigation through side open, curved shaft, 28mm silicone tubing. Irrigation through main hub, aspiration through silicone tubing.


- IA2321 Simcoe Cannula 23G/6#
- IA2323 Simcoe Cannula 22G/7#
- IA2324 Simcoe Cannula 21G/8#



Cannula I/A

IA2322 Reversed Simcoe 23G/6#

Irrigation through silicone tubing, aspiration through main hub.



Cannula I/A 12 O'clock

- IA2325 Vukich 12 O'clock Left I/A Cannula 22G/7#
Angled 90 degrees left 2mm flat tip
- IA2326 Vukich 12 O'clock Right I/A Cannula 22G/7#
Angled 90 degrees right 2mm flat tip



Calipers

Castroviejo Calipers

Measures from 0mm to 20mm in .5mm increments with scale reading on one side.

CA2401 Straight tips

CA2402 Curved tips



Braunstein Fixed Caliper

Fine tips mark sclerotomy sites for vitreous surgery

CA2403 3.5/4.0mm



Braunstein Fixed Caliper

Improved fine tips mark sclerotomy sites for vitreous surgery

CA2409 3.5/4.0mm



Castroviejo Calipers Stainless Steel

Measures from 0mm to 20mm in .5mm increments with scale reading on one side.

CA2411 Straight tips

CA2410 Curved tips



Maloney Astigmatism Keratometer

CA2406 Intraoperative Keratometer, for qualitative measurement of astigmatism.



Punch

Luntz-Dodick Trabeculectomy Punch - Squeeze Handle

- PN2501 0.75mm diameter head
- PN2502 1.0mm diameter head
- PN2503 1.5mm diameter head

Triangle Handle


- PN2504 0.75mm diameter head
- PN2505 1.0mm diameter head
- PN2506 1.5mm diameter head



Trabeculotome Set

Trabeculotome Set left and right

- TB2601 Angled shaft, curved 8mm long pointed with 3mm spread tips (Left)
- TB2602 Angled shaft, curved 8mm long pointed with 3mm spread tips (Right)



Lacrimal Instruments

Lacrimal Dilator

- LD2711 Single-ended, heavy taper
- LD2712 Single-ended, medium taper
- LD2713 Single-ended, fine taper
- LD2714 Castroviejo double-ended, fine/medium taper



Bowman Lacrimal Probe

- LP2721 Bowman Lacrimal Probe, size 0000/000 (0.5/0.6mm)
- LP2722 Bowman Lacrimal Probe, size 00/0 (0.7/0.8mm)
- LP2723 Bowman Lacrimal Probe, size 1/2 (0.9/1.1mm)
- LP2724 Bowman Lacrimal Probe, size 2/3 (1.1/1.3mm)



Lacrimonal Instruments

Pigtail Lacrimonal Probe, *Stainless Steel*

Pigtail curved probes with suture holes just back of blunt ends

PT2715	6mm Ring Diameter
PT2716	7mm Ring Diameter
PT2717	8mm Ring Diameter
PT2718	9mm Ring Diameter



Cannulas


Posterior Capsule Polisher

CA2100	10mm flattened angled and sandblasted tip 18G
CA2101	3mm flattened angled and sandblasted tip 18G



Posterior Capsule Polisher

CA2102	1mm smooth bulbous sandblasted tip 23G
--------	--



Kratz Capsule Scraper

CA2103	1mm angled sandblasted tip 23G
--------	--------------------------------



Irrigating Vectis

CA2104	Thin wall shaft, 4×6mm loop 18G
--------	---------------------------------



Irrigating Vectis

CA2105	Thin wall shaft, 3×8mm loop 18G
--------	---------------------------------



Irrigating Vectis

CA2106	Thin wall shaft, 5×7mm loop 18G
--------	---------------------------------



Cannulas

Irrigating Vectis

CA2107 Thin wall shaft, 3×12mm loop 18G



Knolle Anterior Chamber Irrigating Cannula

CA2126 23G, angled, blunt tip, 2mm from tip to bend, 19mm long excluding hub

CA2127 25G, angled, blunt tip, 2mm from tip to bend, 19mm long excluding hub

Lacrimal Cannula

CA2108 Angled 45° shaft, 15mm from tip to bend
blunt tip, 23G



Lacrimal Cannula

CA2109 Blunt tip, 25G



Lacrimal Cannula

CA2110 45° curved blunt tip, 23G



1×5 Lacrimal Cannula

CA2111 Overall length 72mm, 25G/23G/22G/21G/20G



Bowman Lacrimal Probe

CA2112 0.6mm straight shaft inside, 20G



Cannulas

Lacrimal Cannula

CA2113 10mm with reinforced shaft, 23G



Lacrimal Cannula

CA2114 10mm with reinforced shaft, 23G



West Lacrimal Cannula

CA2115 Gently curved shaft, blunt tip with port, 23G



Nichamin Hydrodissection Cannula

CA2116 Unique flattened tip, 25G



Kellan Hydrodissection Cannula

CA2117 Beveled tip, gently curved, 25G



Retrobulbar Needle

CA2118 Beveled tip 5# 25G

CA2119 Beveled tip 6# 23G

CA2120 Beveled tip 7# 22G

CA2121 Beveled tip 9# 20G



Retrobulbar Needle

28mm from tip to banking pin, with Φ 3mm banking pin

CA2122 Beveled tip 5# 25G

CA2123 Beveled tip 6# 23G

CA2124 Beveled tip 7# 22G

CA2125 Beveled tip 9# 20G



Corneal Ring Set



Optical Zone Marker

KR8301 5.0mm ring
with cross in the center.



Mendez Degree Gauge

MR2801 12mm I.D. 19mm O.D.



Suarez Spreader

KR8302 21mm blade length,
Semi-sharp edge for Epikeratophakia,
115mm



Guider

KR8303 115mm



Double-ended Sinskey Hook

KR8304 Double ended, sinskey
hook and reverse sinskey hook,
110mm



Micro Dissector

KR8305 115mm



Ring Pusher

KR8306 angled tip 3mm, from port,
110mm long



Notched Implantation Forceps

KR8307 0.25mm notched, shafts
angled 45 degrees, with 9mm
platform, 108mm long.



Double Notched Fixation Forceps

KR8308 angled double notched tips,
separated 1mm, 100mm long.



Dissectors


KR8309 5mm Dissectors Set

KR8310 6mm Dissectors Set

Mermoud Set

Mermoud Forceps

MM8321



Mermoud Marker

MM8322 Double end marker 4mm & 5mm



Mermoud Spatula

MM8323



Mermoud Scraper

MM8324



Osher Universal Scissors

CJ3631 Blunt tips, thin curved blades, 21mm from pivot to tip, 115mm length



MISCELLANEOUS

Blade Holder

BB2901 total 115mm



Blade Handle

BH2902 For holding micro blades, 95mm



Bard Parker Handle

BB2916 Handle 125mm with ruler



Chalazion Curette

CR2903 Double-ended, 1.5mm/2.0mm diameter, 130mm length

CR2904 Double-ended 2.5mm/3.5mm diameter, 130mm length



Chalazion Curette

CR2971 Single end, 1.6mm diameter, 110mm length

CR2972 Single end, 2.5mm diameter, 110mm length

CR2973 Single end, 3.0mm diameter, 110mm length



Thornton Fixation Ring

FR2905 Swivel head, 13mm diameter, 9mm cut out

FR2906 Swivel head, 16mm diameter, 14mm cut out



MISCELLANEOUS

Hofmann Thornton Globe Fixation Ring

Improved mechanical structure, better stability and controllability.

FR2965 Swivel head, 13mm diameter full ring

FR2966 Swivel head, 16mm diameter full ring



Hofmann Thornton Globe Fixation Ring

Improved mechanical structure, better stability and controllability.

FR2967 Swivel head, 13mm diameter 3/4 open ring

FR2968 Swivel head, 16mm diameter 3/4 open ring



Kansas Lens Loop

LL2907 2.5*11mm Loop, serrated, for small incision nucleus fragment removal.



Kansas Lens Loop

LL2960 5*7mm Loop



Kansas Lens Loop

LL2961 3*8mm Loop



Kansas Lens Loop

LL2962 4*6mm Loop



MISCELLANEOUS

Wells Enucleation Spoon

ES2908 21mm wide, curved spoon with 5*10mm notch, 147mm long



Cautery Ball

CB2909 Cautery Ball



Anterior Chamber Maintainer

AM2917 Front port O.D. 0.9mm, I.D. 0.6mm
beveled 18G tip, tip length 5.5mm



Anterior Chamber Maintainer

AM2918 3 side ports, O.D. 0.9mm, I.D. 0.6mm,
round 18G tip, tip length 7mm



Anterior Chamber Maintainer

AM2919 Front port, beveled 18G tip,
tip length 3.5mm



Whitman Shepherd Double-Ended Capsule Polisher

CP2911 1.75mm ring diameter, semi-sharp top and bottom edges,
round handle, 140mm long



Berke-Jaeger Lid Plate

LP2915 Curved over length and width, 20mm
and 23mm wide blades, dull finish, 109mm long



MISCELLANEOUS

Hartmann Hemostatic Mosquito Forceps

Ring handle with ratchet lock, 95mm long

HS2921 Straight 20mm serrated jaws

HS2922 Curved 20mm long serrated jaws



Halsted Mosquito Forceps

HS2923 Curved, 32mm from pivot to tip, 125mm long


HS2925 Straight, 32mm from pivot to tip, 125mm long



Halsted Mosquito Forceps Stainless Steel

HS2935 Curved, 32mm from pivot to tip, 128mm long


HS2936 Straight, 32mm from pivot to tip, 128mm long



Halsted Mosquito Forceps Stainless Steel

HS2937 Curved, 42mm from pivot to tip, 160mm long

HS2938 Straight, 42mm from pivot to tip, 160mm long



Schocket Scleral Depressor

Rounded end designed for depressing sclera and curved marking end for reaching behind the globe.

SD2924 2.5mm teardrop shape



Schocket Scleral Depressor

Rounded end designed for depressing sclera and curved marking end for reaching behind the globe.

SD2963 2.5mm teardrop shape



MISCELLANEOUS

Schepens Scleral Depressor

SD2964 17mm diameter finger sleeve, overall length 48mm



Freer Elevator

HS2926 4mm width curved elevator.




Urrets Zavalía Sclera Marker

HS2927 Press and mark the surface of sclera.



Silicone Bulb with cannulas

SB2928 Silicone Bulb with cannulas



Injectors

IOL Injector

IJ2930 Monarch III style Type D



IOL Injector

IJ2931 Monarch III style Type D, plunger mechanism



IOL Injector

IJ2932 Monarch III style Type A B & C, plunger mechanism



Corneal Trephines

Corneal Trephine

17 sizes available 5.0mm through 9.0mm in 0.25mm increments, dull finish

TP2940 Φ 5.0mm	TP 2941 Φ 5.25mm	TP2942 Φ 5.5mm
TP2943 Φ 5.75mm	TP2944 Φ 6.0mm	TP2945 Φ 6.25mm
TP2946 Φ 6.5mm	TP2947 Φ 6.75mm	TP2948 Φ 7.0mm
TP2949 Φ 7.25mm	TP2950 Φ 7.5mm	TP2951 Φ 7.75mm
TP2952 Φ 8.0mm	TP2953 Φ 8.25mm	TP2954 Φ 8.5mm
TP2955 Φ 8.75mm	TP2956 Φ 9.0mm	



Barron Corneal Donor Button Punch

TP2957 For corneal punches 18mm diameter and 14mm height



Sapphire Knives

Phaco Spear

- KN8111 Spear, 60°, 2.2mm, angled
- KN8112 Spear, 60°, 2.4mm, angled
- KN8113 Spear, 60°, 2.6mm, angled
- KN8114 Spear, 60°, 2.75mm, angled
- KN8115 Spear, 60°, 2.8mm, angled
- KN8116 Spear, 90°, 3.0mm, angled

Customized available



Crescent Tunnel

- KN8121 Crescent, 2.0mm, angled
- KN8122 Crescent, 2.2mm, angled
- KN8123 Crescent, 2.8mm, angled



Lance

- KN8131 Lance, 35°, 1.0mm, straight



Sapphire Knives

Single Edge

KN8146 Single Edge 15°, 1.0mm, front cutting



Clear Corneal

KN8139 Clear Corneal, 2.4mm
KN8140 Clear Corneal, 2.5mm
KN8141 Clear Corneal, 2.65mm
KN8142 Clear Corneal, 2.75mm
KN8143 Clear Corneal, 2.8mm
KN8144 Clear Corneal, 3.0mm



Trifacet

KN8145 Trifacet 1.0mm



Trapezoid Spear


KN8147 Trapezoid, 2.0/2.3mm, angled
KN8148 Trapezoid **Customized available**



Sterilizing Trays


Sterilizing Trays

- TR8201 Sterilizing Trays, 160mm*90mm*22mm (Small)
TR8202 Sterilizing Trays, 220mm*140mm*22mm (Middle)
TR8203 Sterilising Trays, 260mm*160mm*22mm (Large)
TR8205 Sterilizing Trays, 160mm*70mm*22mm (Extra Small)




Sterilizing Trays

- TR8204 Sterilizing trays, 265*162*45mm (Double levels)




Sterilizing Trays *Aluminum alloy*

- TR8206 Sterilizing trays, 260mm*155mm* 32mm (Large)



Sterilizing Trays *Aluminum alloy*

- TR8207 Sterilizing trays, 260*160*50mm (Double levels)



20G/23G Forceps

	VF9103/11505	20G Gripping Forceps	1.25mm jaws Straight Shaft, Sharp tips.
	VF9203/21505	23G Gripping Forceps	1.25mm jaws Straight Shaft, Sharp tips.
	VF9104/11506	20G Gripping Forceps	1.25mm jaws Straight Shaft, Round Blunt tips.
	VF9204/21506	23G Gripping Forceps	1.25mm jaws Straight Shaft, Round Blunt tips.
	VF9105/11507	20G Serrated Forceps	1.25mm jaws Straight Shaft
	VF9205/21507	23G Serrated Forceps	1.25mm jaws Straight Shaft
	VF9215/21507C	23G Serrated Forceps	1.25mm jaws Curved Shaft
	VF9106/11508	20G Toothed Forceps	1 x 2 teeth 1.25mm jaws Straight Shaft
	VF9206/21508	23G Toothed Forceps	1 x 2 teeth 1.25mm jaws Straight Shaft
	VF9107/11509	20G Gripping Forceps	with flat Platform 1.25mm jaws Straight Shaft
	VF9207/21509	23G Gripping Forceps	with flat Platform 1.25mm jaws Straight Shaft
	VF9108/11510	20G Membrane peeler	Straight Shaft, 1.25mm Smooth Jaws.
	VF9208/21510	23G Membrane peeler	Straight Shaft, 1.25mm Smooth Jaws.
	VF9109/11511	20G Membrane peeler	Straight Shaft, 1.25mm Grooved Jaws.
	VF9209/21511	23G Membrane peeler	Straight Shaft, 1.25mm Grooved Jaws.
	VF9110/11528	20G Foreign Body Forceps	Three gripping arms Straight Shaft
	VF9111/11512	20G Fine Gripping Forceps	3mm jaws Straight Shaft
	VF9211/21512	23G Fine Gripping Forceps	3mm jaws Straight Shaft
	VF9112/11516	20G Pick and Peel Forceps	3mm Straight Shaft
	VF9212/21516	23G Pick and Peel Forceps	3mm Straight Shaft
	VF9113/11523	20G Pick and Peel Forceps	2.2mm Straight Shaft
	VF9213/21523	23G Pick and Peel Forceps	2.2mm Straight Shaft
	VF9117	ILM End-Gripping Forceps 20G	
	VF9217	ILM End-Gripping Forceps 23G	

Vitreo-Retinal 20G/23G

20G/23G Scissors



VS9301/13151 20G Straight Scissors 1.25mm blade Straight Shaft, Sharp tips.

VS9401/23151 23G Straight Scissors 1.25mm blade Straight Shaft, Sharp tips.



VS9302/13152 20G Straight Scissors 1.25mm blade Straight Shaft, Blunt tips.

VS9402/23152 23G Straight Scissors 1.25mm blade Straight Shaft, Blunt tips.



VS9303/13153 20G Curved Scissors 1.25mm blade Straight Shaft.

VS9403/23153 23G Curved Scissors 1.25mm blade Straight Shaft.



VS9305/13155 20G Vertical Scissors, vertical 40° Straight Shaft

VS9405/23155 23G Vertical Scissors, vertical 40° Straight Shaft



VS9407/23169 23G Curved Scissors 25mm Curved Long Shaft, 3mm blade.



VS9408/23170 23G Curved Scissors 15mm Curved Short Shaft, 2.5mm blade.



Straight Shaft



Curved Shaft

Please refer to the Previous Pages for Photo Details

Retinal Forceps

VF9103	20G Gripping Forceps, Straight Shaft, Sharp tips.
VF9203	23G Gripping Forceps, Straight Shaft, Sharp tips.
VF9104	20G Gripping Forceps, Straight Shaft, Round tips.
VF9204	23G Gripping Forceps, Straight Shaft, Round tips.
VF9105	20G Serrated Forceps, Straight Shaft.
VF9205	23G Serrated Forceps, Straight Shaft.
VF9215	23G Serrated Forceps, Straight Shaft.
VF9106	20G Toothed Forceps, Straight Shaft.
VF9206	23G Toothed Forceps, Straight Shaft.
VF9107	20G Gripping Forceps, Straight Shaft,
VF9207	23G Gripping Forceps, Straight Shaft.
VF9108	20G Membrane peeler, Straight Shaft, Smooth Jaws.
VF9208	23G Membrane peeler, Straight Shaft , Smooth Jaws.
VF9109	20G Membrane peeler, Straight Shaft, Grooved Jaws.
VF9209	23G Membrane peeler, Straight Shaft, Grooved Jaws.
VF9110	20G Foreign Body Forceps, Straight Shaft.
VF9111	20G Fine Gripping Forceps, Straight Shaft.
VF9211	23G Foreign Body Forceps, Straight Shaft.
VF9112	20G Pick and Peel Forceps, Straight Shaft.
VF9212	23G Pick and Peel Forceps, Straight Shaft.
VF9113	20G Pick and Peel Forceps, Straight Shaft.
VF9213	23G Pick and Peel Forceps, Straight Shaft.
VF9117	20G ILM End-Gripping Forceps
VF9217	23G ILM End-Gripping Forceps

Please refer to the Previous Pages for Photo Details

Retinal Scissors

VS9301	20G Straight Scissors, Sharp tips.
VS9401	23G Straight Scissors, Sharp tips.
VS9302	20G Straight Scissors, Blunt tips.
VS9402	23G Straight Scissors, Blunt tips.
VS9303	20G Curved Scissors, Curved blade.
VS9403	23G Curved Scissors, Curved blade.
VS9305	20G Vertical Scissors, vertical 40°
VS9405	23G Vertical Scissors, vertical 40°
VS9407	23G 25mm Curved shaft Scissors.
VS9408	23G 15mm Curved shaft Scissors.

25G Forceps & Scissors



Serrated Gripping Forceps 25G

VF9501 Blunt-tipped forceps with serrated internal surface. Serrated grasping tips for strong membrane grasping. Use for Grasping of strong adherent PVR membranes.



ILM End-Gripping Forceps 25G

VF9502 ILM end-gripping platforms. Improved visualization for grasping thin membranes; Use for ILM- and ERM peeling Fine membrane peeling.



Asymmetrical End-Gripping Forceps 25G

VF9503 Asymmetrical tip allows excellent visualization at grasping area. Interior end-gripping platforms designed for assured grasping; Versatile tip design targets ILM, fine epiretinal membranes & diabetic PVR membranes.



Eckardt End-Gripping Forceps 25G

VF9504 Wider tip enhances surface grasping reduces tearing effect, tip design allows optimal visualization of surgical field; Use for ILM and ERM peeling General peeling of fine membranes which needs end-gripping action.



Kuhn 45° Power Gripping Forceps 25G

VF9505 This Forceps is designed for peeling proliferative intravitreal, epi and sub-retinal membranes, lens capsule, iris, small IOFB, fibrinous membranes and more. The forceps is also suitable for peeling standard ILM and ERM membranes and attached posterior hyaloids.



Vitreo-Retinal 25G

Shah X-tra Grip Forceps 25G

VF9506 Serrations angled 45° for maximum grasping of tissue. Fine tips allow for ILM lifting; Effective for all types of tissue manipulation, ILM, ERM, PVR and diabetic tissue.



Straight Scissors 25G

VS9507 Sharp-tipped straight cutting scissors; For cutting in the iris plane & peripheral membranes.



Curved Scissors 25G

VS9508 Sharp-tipped horizontal curved cutting scissors; Curvature of 12mm to follow the contour of the retina for improved safety; Use for cutting membranes adherent to the retinal surface & All-purpose horizontal membrane cutting.



Micro Incision Instruments

Micro Capsulorhexis Forceps

VC9114 20G Micro Capsulorhexis Forceps


VC9214 23G Micro Capsulorhexis Forceps




Kawai Capsulorhexis Forceps

Serrations, sharp tip can be used as Cystotome. Serrations extend to end of the jaws for better grasp of the anterior capsule in the presence of viscoelastic. Recessed space in jaws prevents inadvertent iris pickup. Jaws open wider for better visibility.

VC9216 23G Micro Capsulorhexis Forceps



Micro Incision Instruments



Micro Graspers Forceps

This utility forceps that can be used for various purposes in the anterior chamber. It can be used for grasping iris to do PI's, performing goniosynechialysis, grasping pupil margins for movement, and more.

VC9218 23G Straight Shaft Forceps

VC9219 23G Curved Shaft Forceps



Coaxial Holder for IOL removal 20G

VF9115 20G coaxial holder for IOL removal (to be used with IOL removal Micro-Scissors VS9015)



Coaxial Micro-Scissors for IOL removal 19G

VS9015 19G coaxial Micro-Scissors for IOL removal (to be used with holder VF9115)





ACMED is one of the leading manufacturers of ophthalmic instruments in China. The manufacturing of ophthalmic surgical instruments traces back to the 1990s. With our meticulous efforts, ACMED are helping thousands of surgeons to cure patients all over the world by our affordable but reliable reusable and disposable surgical instruments. Together with our distributors and partners all over the world, we are committed to provide the best service and guarantee for our products. ACMED also strongly focuses on the innovation of the surgical area. Every day we are continuously receiving new ideas and suggestions from our clients to make ourselves better. We hope to improve and see a better world together with you!



Website: www.acmedglobal.com

Tel/Fax: +86-25-85996718

Email: sales@acmedglobal.com